

# Český metrologický institut

Okružní 31, 638 00 Brno

tel. +420 545 555 111, fax +420 545 222 728, www.cmi.cz

Supervised by:

Regional branch Prague, Radiová 3, 102 00 Prague, Czech Republic,

Department of gas mixtures and reference materials certification - CORM

phone +420 266 020 166, fax +420 272 704 852

## CERTIFICATE

1014-CM-70006-12

### CERTIFIED REFERENCE MATERIAL CZ 70006

POLYCHLORODIBENZO-p-DIOXINS, POLYCHLORODIBENZO FURANS AND TOXIC POLYCHLORO BIPHENYLS IN SEWAGE SLUDGE OF MIXED ORIGIN

Date of issue:

20. December 2006

page 1 / 4

Valid until:

20. December 2011

Date of recertification:

27. February 2012

Valid until:

20. December 2016

Use:

The material is intended for use in verification of analytical procedures and validation of analytical methods and for ensuring traceability. It is not

intended for use as a calibrant.

Description of the material:

The sample consists of 60±0.5g of dried sewage sludge powder (a particle size of less than 100 µm) in amber glass bottles with a polypropylene insert and screw cap. The total I-TEQ value is approximately 150 ng/kg.

The material is radiationally sterilised by <sup>60</sup>Co at a dose of 25 kGy.

Subcontractor:

ANALYTIKA® spol. s r.o., Ke Klíčovu 2a/816,

190 00 Prague 9 Czech Republic

Responsible person:

Dr. Václav Sychra, Ph.D.

Production, testing and certification of this material is in compliance with the Czech Metrology Institute Methodical procedure No. 114-MP-C101 "Preparation and certification of reference materials", and ISO REMCO Guides 34, 35.

Responsible person - CORM:

Director of Regional branch, Prague:



Ing. Vladimír Peršl

Ing. Jan Beránek This Certificate may be published or distributed only in full.

#### Instructions for use:

For analysis the sample should be taken as it is. Its dry mass should be determined on non-analysed portions oven dried at 105°C till the constant mass. The water content is approximately 6% by mass. The recommended sample intake is 5g. The material in the bottle should be rehomogenized before use by manual or mechanical shaking for 1-2 minutes. The bottle should be open several minutes after shaking to avoid possible emission of fine particles of the sample from the bottle and their sedimentation in the surroundings. Storage of the material at room temperature (15-30°C) in the dark is recommended. Since the material is toxic, handling with it should be allowed only to staff trained for working with toxic materials.

#### Certified values and their uncertainties

Analyte		Mass fraction		Number of accepted
		Certified value (ng/kg) <sup>1)3)</sup>	Uncertainty (ng/kg) <sup>2)</sup>	sets of data
2,3,7,8 TeCDD	(D48)	4.5	0.3	7
1,2,3,7,8 PeCDD	(D54)	2.1	0.3	7
1,2,3,4,7,8 HxCDD	(D66)	2.6	0.5	7
1,2,3,6,7,8 HxCDD	(D67)	5.0	0.9	9
1,2,3,7,8,9 HxCDD	(D70)	3.7	0.6	9
1,2,3,4,6,7,8 HpCDD	(D73)	65	10	7
OCDD	(D75)	519	74	7
2,3,7,8 TeCDF	(F83)	110	17	8
1,2,3,7,8 PeCDF	(F94)	157	21	8
2,3,4,7,8 PeCDF	(F114)	87	11	6
1,2,3,4,7,8 HxCDF	(F118)	376	63	9
1,2,3,6,7,8 HxCDF	(F121)	102	13	8
1,2,3,7,8,9 HxCDF	(F124)	11.0	2.2	7
2,3,4,6,7,8 HxCDF	(F130)	19.8	2.8	6
1,2,3,4,6,7,8 HpCDF	(F131)	256	41	9
1,2,3,4,7,8,9 HpCDF	(F134)	110	17	8
OCDF	(F135)	1590	290	9
PCB 77		2380	370	10
PCB 81		108	16	6
PCB 126		169	32	9
PCB 169		25	4	8
PCB 105		3430	500	10
PCB 114		169	36	6
PCB 118		15800	2300	10
PCB 123		121	30	6
PCB 156		9140	1300	8
PCB 157		802	130	10
PCB 167		4130	670	10
PCB 189		1860	260	8

- 1) This value is the unweighed mean of the means of accepted sets of results
- 2) Uncertainties are combined uncertainties multiplied by a coverage factor k=2
- 3) Corrected for dry mass at 105°C

Český metrologický institut Oblastní inspektorát Praha Radiová 3 102 00 Praha

## Non-certified values of contents of some other analytes present in CZ 70006\*

Metals		Polycyclic Aromatic Hydrocarbons (PAHs)		
Analyte	Mass fraction <sup>a)</sup> mg/kg <sup>c)</sup>	Analyte	Mass fraction <sup>a)</sup> mg/kg	
Ag	6.75	acenaphtene	0.53	
As	17.3	anthracene	0.20	
Be	0.82	benzo(a)anthracene	0.61	
Cd	1.47	benzo(b)fluoranthene	0.33	
Co	16.8	benzo(k)fluoranthene	0.21	
Cr	498	benzo(g,h,i)perylene	0.21	
Cu	318	benzo(a)pyrene	0.23	
Hg	1.48	chrysene	0.53	
Mn	1410	dibenzo(a,h)anthracene	0.04	
Mo	4.15	fluoranthene	3.84	
Ni	41.2	fluorene	0.51	
Pb	174	indeno(1,2,3-c,d)pyrene	0.23	
Se	1.41	naphtalene	1.32	
Sn	15.2	phenanthrene	1.92	
Tl	0.17	pyrene	2.27	
V	39.5	Some other organic analytes		
Zn	782	Analyte	Mass fraction <sup>b)</sup> , μg/kg	
Al	(19800) <sup>d)</sup>	$\sum$ HCH, $\alpha$ , $\beta$ , $\gamma$ , $\delta$	(52) <sup>d)</sup>	
Ca	(98200)	2,4'- DDD	(670)	
Fe	(21500)	4,4'- DDD	(1650)	
K	(5650)	2,4'- DDT	(190)	
Mg	(4940)	4,4'- DDT	(850)	
Na	(5470)	1,2,3 – trichlorobenzene	(280)	
		pentachlorobenzene	(200)	
		hexachlorobenzene	(6000)	

Other PCBs ("markers")		Brominated Flame Retardants (BFRs)	
Analyte	Mass fraction <sup>a)</sup> , μg/kg	Analyte	Mass fraction <sup>a)</sup> , μg/kg
PCB - 28	26.3	PBDE - 28	0.77
PCB - 52	14.6	PBDE - 47	62.2
PCB - 101	49.3	PBDE - 49	2.23
PCB - 138	115	PBDE - 66	1.33
PCB - 153	153	PBDE - 85	3.50
PCB - 170	64.8	PBDE - 99	59.8
PCB - 180	123	PBDE - 100	14.0
		PBDE - 153	4.05
		PBDE - 154	4.60
		PBDE - 183	6.25
		PBDE - 209	55.6
		HBCD	79.2

<sup>\*</sup>Analyte contents were derived (unless otherwise stated) from a minimum of 3 data sets

Český metrologický institut Oblastní inspektorát Praha Radiová 3 102 00 Praha

a) Corrected for dry mass at 105°C

b) Corrected for dry mass at 40°C

c) Aqua regia extraction according to ISO-11885

d) Informative analyses only from 2 labs

#### Characterization:

Characterization was performed in accordance with ISO REMCO Guide 35. The certified and non-certified values of analytes were derived from results of an interlaboratory comparison in which 11 selected laboratories took part. Analytical methods used covered an extraction with one or a mixture of solvents and a clean-up in one or several steps. Separation and quantification was carried out by GC-HRMS using two chromatographic phases of different polarity.

#### Metrological traceability:

BCR 614 (S<sub>0</sub>-S<sub>5</sub>), BCR 365 and ERM AC 820-822 were used for calibration by most of participating laboratories. Each laboratory analysed simultaneously with the samples also the matrix reference material BCR – 677.

#### Participating laboratories

#### Material preparation:

ANALYTIKA® spol. s r.o. Prague, Czech Republic

Sub-suppliers:

Ecolab s.r.o., Radlík, Czech Republic (sourcing, preliminary analyses, packing)

Institue of Nuclear Physics, Řež u Prahy (radiational sterilization) Safina a.s. Vestec u Prahy, Czech Republic (homogenization) ÚNS a.s. Kutná Hora, Czech Republic (milling, sieving)

## Homogeneity and stability study:

Axys-Varilab s.r.o., Vrané nad Vltavou - Skochovice, Czech Republic

#### Certification analyses:

- •ALS Czech Republic s.r.o. (former Ecochem a.s.), Prague, Czech Republic
- · Axys Analytical Services Ltd., Sidney, Canada
- Axys-Varilab s.r.o., Vrané nad Vltavou Skochovice, Czech Republic
- •CARSO, Centre d'Analyse et de Recherche sur de Substances Organiques, Lyon, France
- •Institute of Chemical Technology, Metrological Laboratory, Prague, Czech Republic
- Institute of Health, Ostrava, laboratory Frýdek-Místek, Czech Republic
- •MAS, Műnster Analytical Solutions GmbH, Műnster, Germany
- ·Nablabs Laboratories, Espoo, Finland
- •SAL, Scientific Analysis Laboratories Ltd., Manchester, United Kingdom
- •VITO, Vlaamse Instelling voor Technologisch Onderzoek, Mol, Belgium
- Water Research Institute T. G. M., Prague, Czech Republic

#### Note:

A detailed Certification Report on the sample preparation, the analysis procedures and the treatment of the analytical data is available on Internet (www.analytika.net). A paper copy can be obtained from the producer on explicit request.

Český metrologický institut Oblastní inspektorát Praha Radiová 3 102 00 Praha